

# **LARGE FACILITY ULTRA CLEAN AIR COMPRESSORS**

- **TOTALLY ENCLOSED, OSHA APPROVED, BELT GUARD**
- **ASME AIR RECEIVER**
- **LARGE AIR COOLED AFTER COOLERS**
- **LARGE AIR INLET FILTERS AND SILENCERS**
- **ENERGY EFFICIENT MOTORS**
- **MAGNETIC MOTOR STARTERS WITH OVERLOAD PROTECTION**
- **PRESSURE SWITCH OPERATED AUTOMATIC START AND STOP CONTROLS**
- **DISCHARGE AIR, COALESCING FILTER**
- **REFRIGERATED DRYERS**
- **ALL UNITS COME READY TO INSTALL**



***TW* TECH WEST INC.**

Manufacturers of Dental Vacuum  
and Air Systems

2625 N. Argyle Ave. • Fresno, CA 93727  
(559) 291-1650 • (800) 428-7139 • FAX (559) 348-9677

# LARGE FACILITY COMPRESSORS

## FEATURES

Dryer bypass valves.

Finned cylinder and cylinder heads.

Drive pullys.

Special designed flywheel supplies high volume airflow.

All compressors are wired for low voltage.

Crankcase lubrication level indicator.

All compressor tanks are epoxy coated on the inside to meet ASME standards.

## BENEFITS

Provides uninterrupted service during compressor maintenance.

Provides even heat dissipation of the heads.

Dynamically balanced to maximize smooth running conditions.

Cools cylinders and heads.

Allows use of 24V remote switching.

Easy to see the amount and condition of the oil.

To prevent rust and corrosion damage to the tanks.

# LARGE FACILITY COMPRESSORS

MODEL NUMBER	# OF USERS	NO. OF MOTORS x H.P.	VOLTS	AIR DRYER	TANK SIZE	DB(A)	EQUIP. HEIGHT	EQUIP. DEPTH	EQUIP. LENGTH
ACL10S5	10	1 x 5HP	230/460	YES	80	80	47	34	72
ACL20D5	20	2 x 5HP	230/460	YES	120	84	55	46	80
ACL15S5	15	1 x 7.5HP	230/460	YES	120	84	55	40	77
ACL30D5	30	2 x 7.5HP	230/460	YES	200	86	57	50	80
ACL20S5	20	1 x 10HP	230/460	YES	120	86	55	42	77
ACL40D5	40	2 x 10HP	230/460	YES	240	88	58	46	82
ACL30S5	30	1 x 15HP	230/460	YES	120	86	55	40	77
ACL60D5	60	2 x 15HP	230/460	YES	240	88	60	48	86

Revised 1/08

**TECH WEST LARGE FACILITY COMPRESSORS HAVE  
A FULL 2 YEAR WARRANTY**

Exceeds NFPA 99C Guidelines

Tech West's Large Facility Ultra Clean Air Compressors lead the industry in design and quality providing ultra clean, dry air to the dental office. They are built to meet the needs of the most demanding dental practices. They range from 10 to 60 users in a single or dual configuration. They feature our high performance refrigerent dryers. They have air cooled after coolers with heavy cast iron heads for long lasting performance.

At Tech West our large facility Ultra Clean Air Compressor can be custom designed to fit any need or application. Tech West will provide plumbing and electrical design upon request.

